


sck cen



Virtual Conference on Applied  
Radiation Metrology (vCARM) 2021

## Determination of $^{36}\text{Cl}$ and $^{129}\text{I}$ in nuclear solid decommissioning samples

Our challenges:  $^{36}\text{Cl}$  memory effect and how to deal with it?

Inés Llopart Babot<sup>1,2,3</sup>, M. Vasile<sup>1</sup>, A. Dobney<sup>1</sup>, S. Boden<sup>1</sup>, M. Bruggeman<sup>1</sup>, M. Leermakers<sup>2</sup>, J.

Qiao<sup>3</sup>, P. Warwick<sup>4</sup>

<sup>1</sup>Belgian Nuclear Research Centre, SCK CEN, Mol, Belgium

<sup>2</sup>Vrije Universiteit Brussel, VUB, Brussels, Belgium

<sup>3</sup>Technical University of Denmark, DTU, Lyngby, Denmark

<sup>4</sup>GAU-Radioanalytical, University of Southampton, Southampton, UK



# Contents


1. Framework
2. Aim of the study
3. Experimental set-up
4. How to deal with  $^{36}\text{Cl}$  memory effect?
5. Application of the method
6. Detection limit (DL)
7. Conclusions and further work



# Contents


- 1. Framework**
2. Aim of the study
3. Experimental set-up
4. How to deal with  $^{36}\text{Cl}$  memory effect?
5. Application of the method
6. Detection limit (DL)
7. Conclusions and further work

# Decommissioning phase of a nuclear facility




# Radiological waste characterization

## Measurement of the waste



## Current approach



# $^{36}\text{Cl}$ and $^{129}\text{I}$ in decommissioning scenarios

- Significant in terms of *half-life* and *environmental mobility* for **final waste disposal**

$^{36}\text{Cl}$


- Neutron activation natural  $^{35}\text{Cl}$
- $T_{1/2} = 3.02 \cdot 10^5$  year
- $E_{\text{max}} = 709.6$  keV
- In construction materials (as trace element) and in the primary circuit coolant

$^{129}\text{I}$

- Fission product  $^{235}\text{U}$
- $T_{1/2} = 1.57 \cdot 10^7$  year
- $E_{\text{max}} = 154$  keV
- In spent nuclear fuel and in the primary cooling circuit of a nuclear reactor

$\beta$  emitters

Pure fractions





# Contents

1. Framework
- 2. Aim of the study**
3. Experimental set-up
4. How to deal with  $^{36}\text{Cl}$  memory effect?
5. Application of the method
6. Detection limit (DL)
7. Conclusions and further work

# General and specific aims

- To optimize the quantification of  $^{36}\text{Cl}$  and  $^{129}\text{I}$  in solid samples coming from decommissioning scenarios



## CURRENT ISSUES

- ↑ uncertainty on the inventory of DTM
- ↑ DL
- Turn around time (TAT) ↑
- Chemicals needed ↑



## APPLICATION FOR

- Improving radiological characterization
- Reaching low DL to comply with exemption and clearance principles
- Including of the procedure in routine analysis





# Contents

1. Framework
2. Aim of the study
- 3. Experimental set-up**
4. How to deal with  $^{36}\text{Cl}$  memory effect?
5. Application of the method
6. Detection limit (DL)
7. Conclusions and further work

# Step by step

## 1. Sample combustion



**Pyrolyser**

## 2. Separation



**Cl-resin**

## 3. Measurement



**Liquid Scintillation  
Counting (LSC)**

# 1. Combustion

- Pyrolyser from RADDEC is used (procedure adapted from *Warwick et al. 2010*)



Temperature protocols depending on the matrix

Maximum temperature: 900°C

Flow rate: 200 mL/min

No catalyst/quartz beads

Glass connections and quartz tubes and sample boats

Fully Cl release

Carry gases

Flow rate regulation



Avoid plastic surfaces

Cl and I carriers →  $\eta$  quantification

Running blank samples after active samples

## 2. Radiochemical separation

- Based on Zulauf et al. 2010
- **Cl-resin** cartridges + vacuum box are used



# 3. Sample measurement

## Massic activity

### Liquid Scintillation Counting (LSC)

❖ Wallac Quantulus 1220™ **low-level** LS counter

4.9 mL elution medium



15 mL Optiphase Hisafe III

$$A \left( \frac{Bq}{g} \right) = \frac{CPS_{measured}}{C_{eff} * \eta * m_{sample}}$$

## Chemical recovery quantification ( $\eta$ )

**Ion chromatography (IC)**

Stable Cl

**Inductively Coupled Plasma Mass Spectrometry (ICP-MS)**

Stable I

0,1 mL subsample from elution medium

Calibration with Quench Curves

# Optimization of the procedure

## COMBUSTION

- Temperature protocol
- Holding time at maximum temperature
- Amount of the sample
- **Memory effect**

## SEPARATION

- Compatibility with different media
- **Cleaning steps**

## MEASUREMENT

- Counting efficiency calibration
- Homogeneity and stability



**$^{36}\text{Cl}$  contamination in the blanks**



**$^{129}\text{I}$  no contamination in the blanks**



# Contents

1. Framework
2. Aim of the study
3. Experimental set-up
- 4. How to deal with  $^{36}\text{Cl}$  memory effect?**
5. Application of the method
6. Detection limit (DL)
7. Conclusions and further work

# Cross-contamination / memory effect?

- Cross contamination on the vacuum box
  - Blanks + spiked samples
  - Blanks in a different vacuum box
- Memory effect

**No cross-contamination**

Spiked sample with  $^{36}\text{Cl}$



Blank sample




**$^{36}\text{Cl}$  Memory effect**



1. Quartz wool/beads
2. Moist air
3. Bubblers
4. Sample boats




# Quartz wool / beads?



- ↓ CPM when removing quartz wool and quartz beads
- Not possible to remove completely the effect
- Chemical recoveries not significantly affected

**CPM:** counts per minute


# Moist air?



- No differences in CPM while removing the moist air
- Moist air was not affecting  $^{36}\text{Cl}$  memory effect
- No significant differences on the chemical recovery


**CPM:** counts per minute

# Bubblers?



- ↓ CPM when changing the bubbler
- Not possible to remove completely the effect

# Sample boat?




- ↓ CPM when changing the sample boat
- Not possible to remove completely the effect

**CPM:** counts per minute

# Ageing of the sample boats


(Peng and Redfern 2013)



- Variations of quartz structure while  $\uparrow$  temperature  $\uparrow$  CPM while reusing the sample boat
- Memory effect  $\uparrow$  when  $>$  20 times used
- Chemical recovery affected for the reuse

**CPM:** counts per minute

# $^{36}\text{Cl}$ memory effect



- Extra cleaning steps cannot avoid the memory effect
- Reasonable background level between 8-13 CPM

**Average**

**$6.96 \pm 0.69$  CPM**

**CPM:** counts per minute

## Catalyst used in other set-ups

- In some cases catalyst is needed (depending on the matrix and radionuclides) for **different reactions** of the target elements (*Prabir Basu 3<sup>rd</sup> edition 2018*)

Reduction

Ni

(Pearce and Hill 1986)

Oxidation

CuO

(Herod et al. 2014)  
(Hou 2005)

Pt-alumina

(Croudace et al. 2017)


Acceleration  
breakage bond

Fe<sub>2</sub>O<sub>3</sub>

(Du et al. 2020)

- In some cases catalyst **cannot** be removed → included in the pyrolyser set-up
  - Horizontal Split tube furnaces up to 1100 °C (*ENTECH 1918*)

# Catalyst?



- ↓ CPM when removing the catalyst (CuO)
- Not possible to remove completely the effect

Long-term CI memory effect was considered for AMS measurements (*Pavetich et al. 2013*)

**CPM:** counts per minute



# Approach



$^{36}\text{Cl}$  memory effect ~ **Random** effect



$^{36}\text{Cl}$  memory effect **removal**



$^{36}\text{Cl}$  memory effect **corrected**





# Contents

1. Framework
2. Aim of the study
3. Experimental set-up
4. How to deal with  $^{36}\text{Cl}$  memory effect?
- 5. Application of the method**
6. Detection limit (DL)
7. Conclusions and further work



# Application using blank concrete and graphite

## CONCRETE



• Activity measured • Activity spiked

## GRAPHITE



• Activity measured • Activity spiked



# Contents

1. Framework
2. Aim of the study
3. Experimental set-up
4. How to deal with  $^{36}\text{Cl}$  memory effect?
5. Application of the method
- 6. Detection limit (DL)**
7. Conclusions and further work

# Calculation of the detection limit (DL)

- Evaluated based on measurements of several blanks (non-spiked solid samples)
- Calculated according to ISO 11929

Considering  $^{36}\text{Cl}$  memory effect

Parameters (mean values)	$^{36}\text{Cl}$	$^{129}\text{I}$
<b>DL</b>	25 mBq/g	25 mBq/g
<b>Activity concentration values for exemption</b>	1000 mBq/g	10 mBq/g

(Current Basic Safety Standard Directive (96/29/Euratom))



$^{36}\text{Cl}$
45 mBq/g





# Contents

1. Framework
2. Aim of the study
3. Experimental set-up
4. How to deal with  $^{36}\text{Cl}$  memory effect?
5. Application of the method
6. Detection limit (DL)
- 7. Conclusions and further work**

# Conclusions

- ✓ Optimized procedure can be applied to the analyses of **real samples**
- ✓ Optimized procedure can reach the **required limits** by the legislation for  $^{36}\text{Cl}$
- ✓ Catalyst is not needed in this set-up reducing the issues with  $^{36}\text{Cl}$  memory effect



**Blanks are required** after running an active sample (memory effect evaluation)

## Further work


- Investigation of the application to Plastic Scintillating Microspheres (PSm) for  $^{36}\text{Cl}$  determination
- Improvement of iodine quantification in order to reach the limits required by the legislation of  $^{129}\text{I}$  free release



Interlaboratory comparisons and availability of reference materials for these radionuclides **are needed**



sck cen



Virtual Conference on Applied  
Radiation Metrology (vCARM) 2021

**THANK YOU VERY MUCH FOR  
YOUR ATTENTION!**

For further questions do not hesitate to contact me by email  
[ines.llopart@sckcen.be](mailto:ines.llopart@sckcen.be)

## Copyright © SCK CEN - 2021

All property rights and copyright are reserved.

This presentation contains data, information and formats for dedicated use only and may not be communicated, copied, reproduced, distributed or cited without the explicit written permission of SCK CEN.

If this explicit written permission has been obtained, please reference the author, followed by 'by courtesy of SCK CEN'.

Any infringement to this rule is illegal and entitles to claim damages from the infringer, without prejudice to any other right in case of granting a patent or registration in the field of intellectual property.

### **SCK CEN**

Belgian Nuclear Research Centre  
Studiecentrum voor Kernenergie  
Centre d'Etude de l'Energie Nucléaire

Foundation of Public Utility  
Stichting van Openbaar Nut  
Fondation d'Utilité Publique

### **Registered Office:**

Avenue Herrmann-Debrouxlaan 40 - 1160 BRUSSELS - Belgium

### **Research Centres:**

Boeretang 200 - 2400 MOL - Belgium  
Chemin du Cyclotron 6 - 1348 Ottignies-Louvain-la-Neuve - Belgium